

Table Top Nanoindentation Tester (TTX)



The CSM Table Top Instrument combines all the advantages of the standard Nanoindentation Tester (NHT) into a small and simple-to-use instrument which is ideally suited to routine nanoindentation testing where a full platform system may not be appropriate or necessary.

The Table Top instrument has the advantage of being extremely compact and specifically designed for mechanical hardness testing. It can be configured with different motorized or manual tables in order to match restricted budgets. The Table Top instrument therefore provides a cost-effective solution, whilst maintaining the great technical features of the original CSM nanoindentation technique.

The nanoindentation head offers the same high performance as the standard Nanoindentation Tester (NHT) head. It is also upgradeable with optional technical features like the Sinus mode (Dynamic Mechanical Analysis: storage and loss moduli) and can be configured with a range of displacement stage options and integrated video microscope.

Maximum Load	500 mN
Load Resolution	0.04 μ N
Maximum Depth	200 μ m
Depth Resolution	0.04 nm

The TTX consists of the most sophisticated nanoindentation system where all testing modes of nanoindentation are available. Indeed, it offers all the following testing modes:

- Simple indentation
- Advanced indentation
- Progressive multi cycles
- Continuous Multi Cycles (CMC)
- User-defined sequence
- Sinus mode (optional)

CONTACT US FOR A FREE DEMONSTRATION
info@csm-instruments.com +41 32 557 5600

The Table Top Nanoindentation Tester version includes a control unit, a separate Windows XP-based PC and a dual 17" LCD monitor with an acquisition card.